

EN



ElectraSyn 2.0 ElectraSyn 2.0

Introduction

/// ElectraSyn 2.0

In a unique endeavor, IKA owner René Stiegelmann has partnered up with Professor Phil S. Baran of the worldrenowned Scripps Research Institute in La Jolla, California (USA). Over the past three years, engineers and chemists have worked hand in hand to develop a product combining two divisions of chemistry communities: the Electrochemical and the mainstream Synthetic Organic Chemistry Communities that traditionally don't have much in common.

"The culmination of years of research at the crossroads of engineering and synthesis", as Phil Baran describes the development of a new product that has its developmental origins dating back to an apparatus patented and developed by IKA in the 1920s. ElectraSyn 2.0 combines three products in one and will facilitate the mass adoption of electrochemistry for preparative organic synthesis. Society is always positively impacted by the development of sustainable reaction procedures that generate high yields, require less chemical reagents and thus produce less chemical waste. Synthetic organic electrochemistry is an innately sustainable and environmentally friendly field whose widespread adoption has been primarily limited by gaps in engineering rather than desire or potential.

6

/// PACKAGES

21

/// TECHNICAL DATA

22

/// SCALE UP

/// THREE PRODUCTS IN ONE


24

/// IKA SERVICE

- 1. Moeller, K. D. Synthetic Applications of Anodic Electrochemistry. Tetrahedron 2000, 56, 9527–9554.
- 2. Francke, R.; Little, R. D. Redox Catalysis in Organic Electrosynthesis:
- Basic Principles and Recent Developments. Chem. Soc. Rev. 2014, 43, 2492–2521.

 3. Yoshida, J.; Kataoka, K.; Horcajada, R.; Nagaki, A. Modern Strategies in Electroorganic Synthesis. Chem. Rev. 2008, 108, 2265–2299.

 4. Horn, E. J.; Rosen, B. R.; Baran, P. S. Synthetic Organic Electrochemistry: an Enabling and Innately Sustainable Method. ACS Cent. Sci. 2016, 2, 302–308.



The Beginnings: Apparatus for electrolysis with stirring standmotors and regulating starters for the rotating electrodes (electrochemical method)



What is ElectraSyn 2.0?

/// 3 devices in 1
















ElectraSyn 2.0

/// Who is it for?

Discovery scale synthetic organic chemistry: medicines, materials, agrochemicals, natural products, and more.



ElectraSyn 2.0

/// Ready To Go package

This Ready To Go package contains everything you need for your first electrochemical experiment: a potentiostat, a 10 ml glass vial, a holder for the glass vial, a stirring attachment (aluminum) and a set of electrodes (graphite).

This modular equipment concept is ideal for countless different set-ups.

ElectraSyn 2.0 Ready To Go package > ElectraSyn 2.0 base unit + vial holder + single vial, 10 ml, complete + 2 electrodes (graphite) + stir plate (aluminum) + stir bar Ident. No. 0020008980 € 1.749,00 PATENT PENDING

ElectraSyn 2.0 pro

/// CV package







DOZENS OF AVAILABLE ELECTRODES

PREVIOUSLY ANECDOTAL AND BALKANIZED, ELECTROCHEMICAL EXPERIMENTS CAN NOW BE FULLY STANDARDIZED FOR ROBUST, GLOBAL

REPRODUCIBILITY

Graphite SK-50 electrode, 12 pcs. *Ident. No. 0040002858* € **138,00**

Glassy Carbon electrode, 2 pcs. *Ident. No. 0040002842* € **103,00**

Lead bronze electrode, 12 pcs. *Ident. No. 0020016076* € **235,00**

Lead electrode, 12 pcs. *Ident. No. 0040002843* **€ 259,00**

Tungsten electrode, 12 pcs. *Ident. No. 0040002845* **€ 263,00**

Niobium electrode, 12 pcs. *Ident. No. 0040002846* **€ 73,00**

Copper electrode, 12 pcs. *Ident. No. 0040002847* € **72,00**

Magnesium electrode, 12 pcs. *Ident. No. 0040002848* **€ 253,00**

Titanium (grade 2) electrode, 12 pcs. *Ident. No. 0040002849* € **235,00**

Zinc electrode, 12 pcs. *Ident. No. 0040002850* € **235,00**

Stainless Steel electrode, 12 pcs. *Ident. No. 0040002851*

Platinum Plated electrode, 2 pcs. *Ident. No. 0040002852*

€ 89,00 Gold Plated electrode, 2 pcs.

Ident. No. 0040002853 € 237,00

Silver Plated electrode, 2 pcs. *Ident. No. 0040002854* € **44,00**

Aluminum electrode, 12 pcs. Ident. No. 0040003174

Boron Doped Diamond electrode, 2 pcs. *Ident. No. 0040002856* € **410,00**

Tin electrode, 12 pcs. *Ident. No. 0040002857* **€ 247,00**

Nickel electrode, 12 pcs. *Ident. No. 0040002859* **€ 73,00**

RVC electrode, 12 pcs. Ident. No. 0040002860 On request

Nickel Foam electrode, 12 pcs. Ident. No. 0040002861 On request

Cobalt electrode, 2 pcs. *Ident. No. 0040003385*On request

Reference electrodes *Ident. No. 0040002865*€ 98,00



DEFINING INDUSTRY STANDARD

AID IN COMMUNITY
WIDE REPRODUCIBILITY
OF ELECTROCHEMISTRY

1 ml vial complete coming soon

2 ml vial complete coming soon

5 ml vial complete *Ident. No. 0040003171*€ **104,00**

10 ml vial complete *Ident. No. 0040003170* **€ 97,00**

20 ml vial complete *Ident. No. 0040003168* **€ 98,00**

ANALYTICAL CAPABILITIES





Technical data

Voltage reading accuracy	± 16 mV
Current reading accuracy	± 6,2 μA
Minimum voltage step	10 mV



CYCLIC VOLTAMMETRY (CV)

ANALYTICAL MEASUREMENTS AND PREPARATIVE EXPERI-MENTS ON THE SAME DEVICE

CV DATA THAT IS EXPORTABLE

TO YOUR COMPUTER,
VIEWABLE ON YOUR PHONE,
OR EVEN ELECTRASYN WITH
ITS BEAUTIFUL SCREEN.

ElectraSyn 2.0 ElectraSyn 2.0

MODULARITY



ALL KNOWN MODES OF ELECTROCHEMISTRY ARE USABLE AND WITHIN REACH DIVIDED CELL CHEMISTRY CAN BE DONE

Cell type / condition	divided	undivided
	A	A
	oxidation	oxidation
constant current	A A	3 ® /
	reduction	reduction
		· ·
	oxidation	oxidation
constant potential	▼ °	▼
V	reduction	reduction



USER INTERFACE

Never run an electrochemical reaction before? The "Smart Assist" mode will analyze your reaction before running in order to give you a baseline starting point for electrochemical conditions.



ELECTRASYN IS SMART AND COMMUNICATES WITH YOUR REACTION TO HELP YOU

ADVANCED MODE FOR ELECTROCHEMICAL GURUS

USER-FRIENDLY, INTUITIVE AND AESTHETIC INTERFACE FOR RAPID ADOPTION OF ELECTROCHEMISTRY







CLASSIC IKA PLATE BUILT IN INDUSTRY
LEADING
AND PROVEN STIRRING
TECHNOLOGY

ALNICO MAGNET TECHNOLOGY,
FOR EXCELLENT TEMPERATURE STABILITY AND HIGH RESIDUAL INDUCTION



A DEVICE THAT
IMPROVES
OVER TIME
WITH FIRMWARE
UPDATES VIA USB

WILL STAND THE
TEST OF TIME AND
IMPROVE
WITH IT





INDUSTRY LEADING

AND CHEMICAL RESISTANT OPERATION SURFACE

INTEGRATION

OF HIGHLY SENSITIVE CAPACITIVE-TOUCH-ELE-MENTS — COMFORTABLE AND EASY OPERATION INCREDIBLE DISPLAY

TRANSPARENCY

COMPARED TO OTHER TRADITIONAL DISPLAY MATERIALS HARDENED GLASS ENCLOSED,

FAST RESPONSE DISPLAY

FOR MAXIMUM VISIBILITY
AND CHEMICAL RESISTANCE

INDUSTRIAL, BEAUTIFUL,
HIGH-QUALITY
MATERIAL

DESIGNED IN

CALIFORNIA

AND BUILT IN THE

USA



♦ CONNECTIVITY





(S)

WIFI

*

BLUETOOTH










COMING SOON















EASE OF TRANSPORTATION



Technical Data

/// An overview

Potentiostat	
Nominal voltage (input)	48 VCD
Max. current (input)	1.500 mA
Max. input power	40 W
Voltage output	30 / 10 V
Current output	100 mA
EC-Motor rating output	9 W
Speed range	50 – 400 – 1.500 rpm
Setting accuracy speed	10 rpm
Stirring quantity max. per stirring position (H ₂ O)	100 ml
Stirring bar length	10 mm
Operation elements	Capacitive touch / Turning knob
Speed control	Turning knob
Display	TFT
Analog output	no
RS 232 interface	no
USB interface	yes
IP rating	IP 40
Permissible ambient temperature	+5 - + 40 °C
Permissible relative humidity	80 %
Dimensions incl. single-vial Adapter (W \times H \times D)	130 × 150 × 250 mm
Weight	1,4 kg
Power supply	
Input	100 – 240 VAC 1,5 A 50 – 60 Hz
Output	48 VCD 39,84 W LPS (limited power source)
Protection class	II (double insulated)

Scale up

/// ElectraSyn flow

ElectraSyn flow basic is a system for continuous flow electrosynthesis. The heart of this system is the electrosynthesis flow cell – in short, ElectraSyn flow. It consists of two half cells, each equipped with an electrode. By combining similar and dissimilar half cells / electrodes, ElectraSyn flow provides maximum flexibility for research in the field of electrosynthesis. It also enables the laboratory scale production of a variety of products using electrosynthesis.




ElectraSyn flow basic and ElectraSyn flow eco

Ident. No. 0020014266 | Ident. No. 0020014267

€ 9,200.00 |

€ 6,700.00



PACKAGE INCLUDES

- 1 Nine half cells with accessories.
- 2 One nafion-membrane for cell splitting.
- **3** R 104 stand with cell holder and H44 boss head clamp.
- **4** Power supply adapter with cords.
- **5** Peristaltic pump with tubing.*
- **6** Practical carrying / storage case for small components.

Continuously adjustable power supply		
Voltage	0 - 35 V (± 6 mV)	
Current	0 – 1 A (± 50 μA)	
Mains voltage	100, 115 or 230 V (50 Hz / 60Hz)	
Peristaltic pump		
Flow rate per tube	0,01 - 0,61 ml/min	
Overall flow rate	0,02 – 1,22 ml/min	
Inner tubing diameter	0,25 mm	
Mains voltage	100 – 230 V (50 Hz / 60 Hz)	

ElectraSyn 2.0 ElectraSyn 2.0

Customizing Center

It is important that IKA products perform in real laboratory applications. We have a special program of product solutions that are customized to your individual needs. If you cannot find the right device in our standard product range, please send us the details of the specification you need using the online form. Our team will check the feasibility of the specification and offer you a solution.

Please visit www.ika.com/en to have look at the product modification requests that we have already implemented.

Worldwide service network

/// Direct contacts in your region

Our dedicated team of engineers provides comprehensive technical service on a global level. If you have any questions, please do not hesitate to contact IKA directly. Alternatively, you can get in touch with your dealer.

IKA guarantees that spare parts will be available for 10 years. In the event of any faults with a device, or if you have any technical questions regarding our products, their maintenance or replacement parts, please call us at **00 8000 4524357** (00 8000 IKAHELP) or send an eMail to **service@ika.de**

IKA Application Support

Our Application Center spans 400 m² and is equipped with the most modern facilities for presenting and testing laboratory equipment and processes. The Center brings us even closer to our customers and improves our service. If you are interested, you can use our facilities to test processes that include stirring, shaking, dispersing, grinding, heating, analysis and distillation.

Call us on **00 8000 4522777** (00 8000 IKAAPPS) or send an eMail to applicationsupport@ika.de.



Our Application Center covers 400 m² and offers modern equipment for demonstrating and testing laboratory devices and processes.



Send us your sample. We will run a test with the suitable device within 48 hours.



We would be happy to help you find the **perfect device** for your application.



Interested individuals and customers can **test processes** including stirring, shaking, dispersing, milling, heating, analyzing and distilling.





EN

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany Phone: +49 7633 831-0, Fax: +49 7633 831-98

eMail: sales@ika.de

/// WORLDWIDE

USA

IKA Works, Inc. Phone: +1 910 452-7059 eMail: sales@ika.net

MALAYSIA

IKA Works (Asia) Sdn Bhd Phone: +60 3 6099-5666 eMail: sales.lab@ika.my

JAPAN

IKA Japan K.K. Phone: +81 6 6730 6781 eMail: info_japan@ika.ne.jp

KOREA

IKA Korea Ltd. Phone: +82 2 2136 6800 eMail: info@ika.kr

CHINA

IKA Works Guangzhou Phone: +86 20 8222 6771 eMail: info@ika.cn

INDIA

IKA India Private Limited Phone: +91 80 26253 900 eMail: info@ika.in

BRASIL

IKA do Brasil Phone: +55 19 3772-9600 eMail: sales@ika.net.br

POLAND

IKA Poland Sp. z o.o. Phone: +48 22 201 99 79 eMail: sales.poland@ika.com

ENGLAND

IKA England LTD.
Phone: +44 1865 986 162
eMail: sales.england@ika.com

/// PRE-ORDER NOW

For more information on our products or to order online, please visit: www.ika.com/en

/// SOCIAL MEDIA



IKAworldwide



IKAworldwide /// #lookattheblue



@IKAworldwide

designed by FormulaWorks – made by IKA Works

Subject to technical changes and alteration of prices. Prices valid until December 31, 2017.
Prices exclude taxes and other sales related fees.