

CHROMAFIL® Xtra



Pure Filters

MACHERY-NAGEL

www.mn-net.com



Since 1911



CHROMAFIL® - Disposable filters from MACHERY-NAGEL



Contents

➤ Introduction	3
➤ Technical information	5-9
➤ Specifications	11
➤ Quality control	12
➤ Packages and membranes	13-14
➤ Chemical compatibility	15
➤ Optimal use of CHROMAFIL®	16
➤ How to select the optimal CHROMAFIL®	17
➤ Ordering information	18-29
➤ Technical support and customer service	31



CHROMAFIL®

Ideal for GC, HPLC and UHPLC sample clarification

Introduction

CHROMAFIL® syringe filters are used for filtration of suspended matter from liquid samples (1–100 mL) or gases. With CHROMAFIL®, rapid purification and removal of particles is very simple: just place the filter on the syringe, and you are ready for filtration. Special manipulations are not required. Contamination of sensitive instrumentation by solid impurities can be avoided, thus increasing lifetime of chromatographic columns and equipment. The filter can be used for the sample preparation for HPLC, GC, ICP, AAS, TOC, DOC, IR, NMR, photometry, spectroscopy, . . .

- different membrane types to meet multiple filtration needs
- low content of extractable compounds ensure reliable analyses
- superior chromatography column protection helps extend column life
- fast flow geometry for easy filtration
- low hold-up volume for maximum filtrate recovery
- HPLC certified
- designed to be compatible for use on all common automated filtration systems, e.g. SOTAX dissolution systems


CHROMAFIL® Xtra

labeled for method validation and certification

Xtra imprint for direct identification of the membrane type, diameter and pore size

Xtra low bleeding PP housing

Xtra color-free plain polypropylene







Technical Information

❖ **Low content of extractable substances** due to a high density polypropylene housing combined with ultrapure filtration membranes.

In comparison to filters made of polycarbonate, polyacrylate or polystyrene, all CHROMAFIL® filters are resistant against nearly all organic solvents.

(see list of chemical compatibility on page 15)

HPLC-test

Conditions: 2 mL of the solvent (specified on top of the chromatograms) were applied to the filter;
100 µL of the filtrate were injected to the HPLC.

Eluent A: water

Eluent B: acetonitrile

Gradient: 10 % → 95 % B in 25 minutes


Flow rate: 0.5 mL/min

Sensitivity: -5 to 10 mAU at 254 nm


Column: 125/4 NUCLEODUR® C₁₈ Gravity 5 µm
(REF: 760100.40)

Acetonitrile:

Eluate:



Blanc:







CHROMAFIL® - Disposable filters from MACHEREY-NAGEL

Low Bleeding PP housing

Even a treatment with very aggressive solvents/reagents does not lead to significant blind values or extractables.

For a proof, a filter was rinsed with 0.5 mL MSTFA (N-methyl-N-trimethylsilyl-trifluoroacetamide), a very powerful silylation reagent. The result is shown in the GC/MS chromatogram.





To provide **the lowest content of extractable substances**, the housing of every CHROMAFIL® filter is ultrasonically sealed

- ✦ The filters are **welded, not glued**, because glue may have extractable ingredients
- ✦ The welding leads to a tight connection between the both parts, thus the filter can be used in both directions. No fluid can leak from the filter housing.

The special **thick rim** of the housing is ideal for use in laboratory robots (e. g. SOTAX, Benchmate™, . . .).



For a **safe connection on the „high pressure“ side** every CHROMAFIL® filter provides a **Luer lock on inlet**.



Luer outlet

- ✦ For the 3, 13 and 25 mm diameter filter: standard luer outlet
This luer configuration offers low hold-up volume and easy filtration into autosampler vials and NMR tubes

Filter inlet and filter outlet can be fitted to all CHROMABOND® columns and accessories for selective sample preparation with the aid of a special adaptor.



CHROMAFIL® - Disposable filters from MACHERY-NAGEL



No breakage of the membrane due to a stabilizing “crash” plate

- The sample fluid is deviated in four lanes by the “crash” plate and does not directly hit the membrane. The resulting pressure distribution protects the membrane against breakage.


Optimal flow geometry by star-shaped distribution plate

- The fluid penetrates the membrane on the whole surface, not only on a small area; the **filter will not clog rapidly**, which guarantees in a high flow efficiency.



Different pore sizes for multiple filtration application


Available **pore sizes** 0.2 and 0.45 μm (additional: PET filters with 1.2 μm , glass fiber filters with 1 μm , PES filters with 5 μm). Filters with 0.45 μm pore size remove fine particles which can clog chromatography columns. **0.2 μm pore size filters are recommended for filtration of UHPLC samples.**



Filter Sizes

3, 13 and 25 mm effective membrane diameter. The small diameter filters are especially recommended for very small samples, which require extremely low dead volumes.

Sample volume	Recommended membrane diameter	Dead volume	Filtration area
≤ 1 mL	3 mm	5 µL	0.07 cm ²
1-5 mL	13 mm	30 µL	1.33 cm ²
5-100 mL	25 mm	80 µL	4.91 cm ²





All filters can be **autoclaved** at **121 °C** and **1.1 bar** for 30 min.

25 mm CHROMAFIL® filters are designed to be **100% compatible and reliable for use with the SOTAX AT70 smart** fully automated dissolution testing systems.



CHROMAFIL® - Disposable filters from MACHERY-NAGEL



A specification data sheet

is available for all membranes and filter diameters

Enhanced quality control for better results

MN certifies CHROMAFIL® syringe filters to be low in UV absorbing extractables.

All filters and membrane types have been HPLC tested for compatibility with the most common HPLC solvents (methanol, water, acetonitrile, see test chromatograms).

HPLC-test certificates are available for every membrane type.

Please visit: www.mn-net.com

Example of specification data sheet



CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



Pressure stability of CHROMAFIL® syringe filter housing 12 bar

The „blue“ test: membrane, pressure and filtration batch test with blue colored silica particles in matching particle sizes provides an excellent method to find leaks or membrane deviations.



Package sizes

packs of 100 or 400 (BigBoxes) for 25 mm Ø filter
 packs of 100 for 13 mm Ø filter
 packs of 100 for 3 mm Ø filter
 packs of 50 for sterile filter

Different membrane materials for multiple filtration applications

Depending on your filtration task you can choose filter membranes made from different materials:

Polyester	(PET) with or without glass fiber prefilter*
Regenerated cellulose	(RC) with or without glass fiber prefilter*
Polytetrafluoroethylene	(PTFE)
Hydrophilized Polytetrafluoroethylene	(H-PTFE)
Cellulose mixed esters	(MV)
Cellulose acetate	(CA) sterile and non-sterile
Polyamide / Nylon	(PA)
Polyethersulfone	(PES)
Polyvinylidene difluoride	(PVDF) with or without glass fiber prefilter*
Glass fiber	(GF)

* Filters with (nom. 1 µm) GF prefilter provide a 2-4 times greater throughput than filter without prefilter for extremely viscous and most difficult-to-filter samples



CHROMAFIL® - Disposable filters from MACHERY-NAGEL



Combi syringe filters with a coarse glass fiber (GF) prefilter and a small-pore membrane as main filter

User benefits:

- for solutions with a high load of particulate matter: lower back pressure, easy filtration
- for high yields of filtrate: more mL of pure filtrate per filter

The technology:

The glass fiber membrane (1 μm) removes coarse particles, before they can block the fine main membrane. This results in a better filtration efficiency, particularly for highly contaminated samples.

Housing: solvent-resistant, ultra low bleed polypropylene

Inlet: Luer lock

Outlet: Luer

Pore diameter: 1.0/0.20 μm or 1.0/0.45 μm

Filter diameter: 25 mm

Void volume: < 80 μL

Packing unit: 100 filters / BigBoxes with 400 filters

Available membranes with GF-prefilter:

Polyester (PET)

Regenerated cellulose (RC)

Polyvinylidene Difluoride (PVDF)

CHROMAFIL® - Disposable filters from MACHERY-NAGEL



Solvent	Material										
	MV	CA	RC	PA	PTFE	H-PTFE	PVDF	PES	PET	GF	PP
Acetaldehyde	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Acetic acid, 100 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Acetone	⊖	⊖	⊕	⊕	⊕	⊕	⊖	⊖	⊕	⊕	⊕
Acetonitrile	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Ammonia, 25 %	⊖	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Benzene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
n-Butanol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Cyclohexane	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Dichloromethane	⊕	⊖	⊖	⊖	⊕	⊕	⊕	⊖	⊕	⊕	⊖
Diethyl ether	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Dimethylformamide	⊖	⊖	⊕	⊕	⊕	⊕	⊖	⊖	⊕	⊕	⊕
1,4-Dioxane	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊖	⊕	⊕	⊕
Ethanol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Ethyl acetate	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Ethylene glycol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Formic acid, 100 %	⊕	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Hydrochloric acid, 30 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Methanol	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Nitric acid, 65 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊖
Oxalic acid, 10 % aqueous	⊕	⊕	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Petroleum ether	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Phosphoric acid, 80 %	⊖	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Potassium hydroxide, 1 mol/L	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
2-Propanol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Sodium hydroxide, 1 mol/L	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Tetrachloromethane	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Tetrahydrofuran	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊖	⊕	⊕	⊕
Toluene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Trichloroethylene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Trichloromethane (Chloroform)	⊕	⊖	⊕	⊖	⊕	⊕	⊕	⊖	⊕	⊕	⊖
Urea	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Water	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Xylene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕

The table lists the chemical compatibility of our CHROMAFIL® materials. The chemical compatibility depends on several parameters such as time, pressure, temperature and concentration. In most cases, CHROMAFIL® filters will have only short contact with a solvent. In these cases they may be used despite of limited compatibility. For example, a PTFE filter with PP housing does not release any UV-detectable substances during filtration of 5 mL THF, although PP shows only limited resistance towards THF.

Data not guaranteed.

⊕ resistant, ⊖ not resistant, ⊕ limited resistance

MV = cellulose mixed esters

CA = cellulose acetate

RC = regenerated cellulose

PA = polyamide (Nylon)

PTFE = polytetrafluoroethylene

H-PTFE = hydrophilized polytetrafluoroethylene

PVDF = polyvinylidene difluoride

PES = polyethersulfone

PET = polyester

GF = glass fiber

PP = polypropylene (housing material)



CHROMAFIL® - Disposable filters from MACHERY-NAGEL

Optimal use of CHROMAFIL® syringe filter

For achieving the full benefits of filtration we recommend the following instructions.



Draw up the sample into the syringe. Then draw approximately 1 mL of air into the syringe. The air helps to minimize the remaining fluid in the filter.



Plug the CHROMAFIL® syringe filter onto the syringe with the luer connection. Ensure a tight connection by gently turning.



Start with gentle pressure to filter your sample into a vial*. This helps to assure maximum throughput.



Tips / additional information

We recommend either discarding the first 1 mL or rinsing the filter unit with 1 mL of primary solvent before sample filtration.

In order to avoid breakage of the membrane only syringes with volumes of 10 mL or higher should be used.



Do not reuse syringe filters
Do not use at temperatures above 55 °C (131 °F)

Warning: CHROMAFIL® syringe filters are intended for laboratory use only. Do not use CHROMAFIL® syringe filters for direct patient care applications.

* MACHERY-NAGEL offers a wide range of vials and caps. More information at www.mn-net.com/vials



How to select the optimal CHROMAFIL® syringe filter

1. Filter size		2. Pore size of filter membrane	3. Membrane type
sample volume	filter size	pore size	properties of sample
 5-100 mL	 25 mm	 0.45 µm	aqueous, polar hydrophilic low particle-load high particle-load, prefiltration required mid-polar e.g. HPLC eluents proteins protein needed remove protein strong acids and bases organic, nonpolar, hydrophobic low particle-load high particle-load, prefiltration required
 1-5 mL	 13 mm	 0.20 µm	PET H-PTEF MV RC GF/PET GF/RC GF/PVDF PET PA RC CA PVDF PES GF GF/PET GF/PVDF H-PTEF PTFE PTFE PET GF/PET GF/PVDF
 < 1 mL	 3 mm		



CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



Polyester (PET)

- ✦ hydrophilic multipurpose membrane
- ✦ for polar as well as nonpolar solvents
the HPLC filter, especially suited for mixtures of water and organic solvents for TOC/DOC determination, not cytotoxic, does not inhibit the growth of microorganisms and higher cells
- ✦ PET filters with integrated glass fiber prefilter (**GF/PET**) are recommended for solutions with a high load of particulate matter or for highly viscous solutions

Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters / pack	REF
	Xtra PET-20/25	0.20	25	labeled		100	729221	400	729221.400
	Xtra PET-45/25	0.45	25	labeled		100	729220	400	729220.400
	Xtra PET-120/25	1.2	25	labeled		100	729229	400	729229.400
	Xtra PET-20/13	0.20	13	labeled		100	729222		
	Xtra PET-45/13	0.45	13	labeled		100	729223		
Combi Filters									
	GF/PET-20/25	1.0/0.20	25	blue	orange	100	729032	400	729032.400
	GF/PET-45/25	1.0/0.45	25	black	orange	100	729033	400	729033.400



Regenerated Cellulose (RC)

- ✧ hydrophilic membrane with very low adsorption
- ✧ for aqueous and organic / aqueous liquids, i.e. polar and medium polar sample solutions
- ✧ binding capacity for proteins 84 µg per 25 mm filter
- ✧ RC filters with integrated glass fiber prefilter* (**GF / RC**) are recommended for solutions with a high load of particulate matter or for highly viscous solutions

Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters/pack	REF
	Xtra RC-20/25	0.20	25	labeled		100	729230	400	729230.400
	Xtra RC-45/25	0.45	25	labeled		100	729231	400	729231.400
	Xtra RC-20/13	0.20	13	labeled		100	729236		
	Xtra RC-45/13	0.45	13	labeled		100	729237		
Combi filters									
	GF/RC-20/25	1.0/0.20	25	blue	blue	100	729050	400	729050.400
	GF/RC-45/25	1.0/0.45	25	black	blue	100	729051	400	729051.400

* glass fiber exhibits a high protein-binding capacity



CHROMAFIL® - Disposable filters from MACHERY-NAGEL



Polytetrafluoroethylene (PTFE)

- ✦ hydrophobic membrane
- ✦ for nonpolar liquids and gases
- ✦ very resistant towards all kinds of solvents as well as acids and bases
flushing with alcohol, followed by water, makes the originally hydrophobic membrane more hydrophilic

Ordering information · CHROMAFIL®

	Type	Pore size [μm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra PTFE-20/25	0.20	25	labeled	100	729207	400	729207.400
	Xtra PTFE-45/25	0.45	25	labeled	100	729205	400	729205.400
	Xtra PTFE-20/13	0.20	13	labeled	100	729208		
	Xtra PTFE-45/13	0.45	13	labeled	100	729209		
	O-20/3	0.20	3		100	729014		
	O-45/3	0.45	3		100	729015		



Hydrophilized polytetrafluoroethylene (H-PTFE)

- hydrophobic membrane with additional hydrophilic properties
- for polar and nonpolar sample solutions
- resistant towards all kinds of solvents as well as acids and bases

Ordering information · CHROMAFIL®

	Type	Pore size [μm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra H-PTFE-20/25	0.20	25	labeled	100	729245	400	729245.400
	Xtra H-PTFE-45/25	0.45	25	labeled	100	729246	400	729246.400
	Xtra H-PTFE-20/13	0.20	13	labeled	100	729256		
	Xtra H-PTFE-20/13	0.45	13	labeled	100	729257		



CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



Cellulose Mixed Ester (MV)

- hydrophilic membrane with very low adsorption
- for aqueous or polar solutions

Ordering information · CHROMAFIL®

	Type	Pore size [μm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra MV-20/25	0.20	25	labeled	100	729206	400	729206.400
	Xtra MV-45/25	0.45	25	labeled	100	729204	400	729204.400



Cellulose Acetate (CA)

- ✦ hydrophilic membrane
- ✦ for filtration of water-soluble oligomers and polymers, especially suited for biological macromolecules
- ✦ very high stability in aqueous solutions
- ✦ binding capacity for proteins 21 µg per 25 mm filter
- ✦ also available in a sterile package (S) for filtration under sterile conditions (each filter individually sealed)

Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters / pack	REF
	Xtra CA-20/25	0.20	25	labeled		100	729226	400	729226.400
	Xtra CA-45/25	0.45	25	labeled		100	729227	400	729227.400
	Xtra CA-20/13	0.20	13	labeled		100	729254		
	Xtra CA-45/13	0.45	13	labeled		100	729255		
Sterile filters									
	CA-20/25 (S)	0.20	25	yellow	red	50	729024		
	CA-45/25 (S)	0.45	25	colorless	red	50	729025		



CHROMAFIL® - Disposable filters from MACHERY-NAGEL



Polyamide (PA) = Nylon

- moderately hydrophilic membrane
- for aqueous and organic / aqueous medium polar liquids

Ordering information · CHROMAFIL®

	Type	Pore size [μm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra PA-20/25	0.20	25	labeled	100	729212	400	729212.400
	Xtra PA-45/25	0.45	25	labeled	100	729213	400	729213.400
	Xtra PA-20/13	0.20	13	labeled	100	729248		
	Xtra PA-45/13	0.45	13	labeled	100	729249		
	AO-20/3	0.20	3		100	729010		
	AO-45/3	0.45	3		100	729011		



Polyethersulfone (PES)

- ✦ hydrophilic membrane
- ✦ for aqueous and slightly organic liquids with higher flow rates
- ✦ very low adsorption for pharmaceuticals and proteins
- ✦ good stability against organic acids and bases
- ✦ binding capacity for proteins 29 µg per 25 mm filter

Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra PES-20/25	0.20	25	labeled	100	729240	400	729240.400
	Xtra PES-45/25	0.45	25	labeled	100	729241	400	729241.400
	Xtra PES-500/25	5.0	25	labeled	100	729242	400	729242.400



CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



Polyvinylidene Difluoride (PVDF)

- ✦ hydrophilic membrane
- ✦ for aqueous solutions, water-soluble oligomers and polymers like proteins
- ✦ low binding capacity for proteins 20 µg per 25 mm filter
- ✦ PVDF filters with integrated glass fiber prefilter* (**GF / PVDF**) are recommended for filtration of biological samples with high particle loads.

Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters / pack	REF
	Xtra PVDF-20/25	0.20	25	labeled		100	729218	400	729218.400
	Xtra PVDF-45/25	0.45	25	labeled		100	729219	400	729219.400
	Xtra PVDF-20/13	0.20	13	labeled		100	729243		
	Xtra PVDF-45/13	0.45	13	labeled		100	729244		
Combi filters									
	GF/PVDF-45/25	1.0/0.45	25	black	white	100	729039	400	729039.400



* glass fiber exhibits a high protein-binding capacity



Glass Fiber (GF)

- inert filter, nominal pore size 1 µm, allows higher flow rates than small pore filters
- for solutions with high loads of particulate matter or for highly viscous solutions (e.g. soil samples, fermentation broths)
- as prefilters for other CHROMAFIL® filters, they prevent clogging of the membrane

Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra GF-100/25	nom. 1.0	25	labeled	100	729228	400	729228.400
	GF-100/13	nom. 1.0	13	labeled	100	729234		




CHROMAFIL® - Disposable filters from MACHERY-NAGEL




Filtration cartridges

- ✦ Filtration cartridges for sample clarification under vacuum (e.g., using the CHROMABOND® vacuum manifold or SPE automation systems like Gilson Aspec™, Rapidtrace) or by gravity flow.
- ✦ Cartridge sizes 3 mL and 6 mL
- ✦ Different membranes (PET, RC, PTFE, PVDF, GF) and pore sizes (0.2, 0.45 and 1.0 µm). The membrane materials correspond to the respective CHROMAFIL® syringe filters.

Ordering information · CHROMAFIL® filtration cartridges


	Type	Pore size [µm]	Pack of	Column volume REF	
				3 mL	6 mL
	PET (polyester)	0.20	100	730578.320	730578.620
	PET (polyester)	0.45	100	730578.345	730578.645
	RC (regenerated cellulose)	0.20	100	730068.320	730068.620
	RC (regenerated cellulose)	0.45	100	730068.345	730068.645
	PTFE (polytetrafluoroethylene)	0.20	100	730570.320	730570.620
	PTFE (polytetrafluoroethylene)	0.45	100	730570.345	730570.645
	PVDF (polyvinylidene difluoride)	0.20	100	730579.320	730579.620
	PVDF (polyvinylidene difluoride)	0.45	100	730579.345	730579.645
	GF (glass fiber)	nom. 1.0	100	730517.3100	730517.6100



MULTI 96 filter plates

- 96-well polypropylene plates for simultaneous filtration of 96 samples
- Advantages of this high-throughput system:
Economical by saving time and solvent
Use of multi-channel pipettors facilitates liquid transfer steps
Readily adaptable to all common automated/robotic handling systems
Minimized dead volume ($\leq 40 \mu\text{L}$)
- Membrane materials correspond to the respective CHROMAFIL® syringe filters

Ordering information · CHROMAFIL® MULTI 96 filter plates

	Description	Pack of	REF
	Filter plates with cellulose mixed ester filter elements ($0.20 \mu\text{m}$)	1	738770.M
	Filter plates with cellulose mixed ester filter elements ($0.45 \mu\text{m}$)	1	738771.M
	Filter plates with RC filter elements (regenerated cellulose, $0.20 \mu\text{m}$)	1	738656.M
	Filter plates with RC filter elements (regenerated cellulose, $0.45 \mu\text{m}$)	1	738657.M
	Filter plates with PTFE filter elements ($0.20 \mu\text{m}$)	1	738660.M
	Filter plates with PTFE filter elements ($0.45 \mu\text{m}$)	1	738661.M
	Filter plates with PTFE filter elements ($1.0 \mu\text{m}$)	1	738662.M
	Filter plates with PTFE filter elements ($3.0 \mu\text{m}$)	1	738663.M
	Filter plates with PE filter elements ($20 \mu\text{m}$)	1	738655.M
	Filter plates with PE filter elements ($50 \mu\text{m}$)	1	738659.M
	Filter plates with glass fiber filter elements (nominal $1 \mu\text{m}$)	1	738655.2M
	Filter plates with glass fiber filter elements (nominal $3 \mu\text{m}$)	1	738658.M





Technical Support and Customer Service

... we Meet your Needs

If you have any questions concerning CHROMAFIL® filters or our chromatography program, or if you are looking for solutions to a special application, please feel free to contact us:

www.mn-net.com · info@mn-net.com

Our website offers **more than 3000 applications** which might already solve your analytical questions.

Please visit: www.mn-net.com/apps



**Germany
and international**

phone +49 24 21 969-0
toll-free 0800 2616 000
fax +49 24 21 969-199 or -198
e-mail sales-de@mn-net.com



USA

phone +1 484 821 0984
toll-free 888-321-6224 (MACH)
fax +1 484 821 1272
e-mail sales-us@mn-net.com



France

phone +33 388 68 22 68
fax +33 388 51 76 88
e-mail sales-fr@mn-net.com



Switzerland

phone +41 62 388 55 00
fax +41 62 388 55 05
e-mail sales-ch@mn-net.com

**“FilterFinder” online:
THE cross reference for syringe filter
www.mn-net.com/filterfinder**





HPLC



GC



TLC




SPE and Flash



Syringe filters



Vials and caps



... we Meet your Needs



www.mn-net.com

MACHERY-NAGEL



MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

Germany and international:

Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

Switzerland:

MACHERY-NAGEL AG
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

France:

MACHERY-NAGEL EURL
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

USA:

MACHERY-NAGEL Inc.
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com



Since 1911